



Target group

Presbyopes who spend most of the day indoors, working on a laptop, reading various documents, or talking to clients. Age 45+.

Note: The use of this type of lens in traffic is strictly prohibited.

Technology

OPTOTECH OFFICE 16

Method of production

The internal manufacturing of the digressive surface.

Diopters ordering

1/4

Progressive length

18 mm

Minimum fitting height

20 mm

Method of ordering

Diopter for proximity


Digression selection based on the table (below) or as desired

Use pupillary distance for distance, preferably for each eye individually (IPD)

Optimization

Special thinning for the lens in the plus range based on the shape of the frame and the parameter for fitting the lens to the frame. If the lens is ordered without the Optimization option, the topography of the lens will be decentered nasally in relation to the geometric center by 2.5 mm.

User ratings




Markings on the lenses

engravings
 markers

DIGRESSION SELECTION

Addition	Distance to 1 m	Distance to 2 m	Distance to 3 m
1,00	-	-	0,75
1,25	-	0,75	1,00
1,50	0,75	1,00	1,25
1,75	1,00	1,25	1,50
2,00	1,25	1,50	1,75
2,25	1,50	1,75	2,00
2,50	1,75	2,00	2,25
2,75	2,00	2,25	-
3,00	2,25	-	-



MATERIALS

1.50 CR-39	1
1.53 Trivex	2
1.56 Mid. index	3
1.59 Poly	4
1.60 High Index	5
1.60 Tribid	6
1.67 High Index	7
1.74 High Index	8

OPTIONS

None	W	Transitions GEN 8 Sapphire	F	Nupolar Gray	S
Remove 400	H	Transitions GEN 8 Amethyst	M	Nupolar Gray Light	L
Remove 420	R	Transitions GEN 8 Amber	A	Nupolar Green	E
Remove 420 Photomatic Brown	O	Transitions GEN 8 Emerald	E	Transitions XTRActive Polarised Gray	X
Remove 420 Photomatic Gray	U	Transitions XTRActive Brown	C	Photomatic Brown	B
Transitions GEN 8 Brown	Z	Transitions XTRActive Gray	L	Photomatic Gray	I
Transitions GEN 8 Gray	N	Transitions XTRActive Green	K	DriveWear	D
Transitions GEN 8 Green	V	Nupolar Brown	B		