

MONOFOKAL

HARMONY



Target group

Monofocal lens with a slightly reinforced insert in the lower part of the lens. It is recommended for the category of people who feel strain, tension in the eyes, and headaches due to prolonged work with digital devices.

Design 0.50 – age 20+
Design 0.75 – age 30+

Technology

OPTOTECH Mild Add

Method of production

Internal manufacturing of monofocal surface.

Diopters ordering

1/4

Minimum fitting height

18 mm

Optimization

Special thinning for the lens in the plus range based on the shape of the frame and the parameter for fitting the lens to the frame. If the lens is ordered without the Optimization option, the topography of the lens will be decentered nasally in relation to the geometric center by 2.5 mm.

Target group

Vision – distance



Vision – near vision field



Spontaneous adaptation



Vision – medium distance



Comfort




Use of digital devices



Markings on the lenses

█ engravings
█ markers



MATERIALS

1.50 CR-39	1
1.53 Trivex	2
1.56 Mid. index	3
1.59 Poly	4
1.60 High Index	5
1.60 Tribid	6
1.67 High Index	7
1.74 High Index	8

OPTIONS

None	W	Transitions GEN 8 Sapphire	F	Nupolar Gray	S
Remove 400	H	Transitions GEN 8 Amethyst	M	Nupolar Gray Light	L
Remove 420	R	Transitions GEN 8 Amber	A	Nupolar Green	E
Remove 420 Photomatic Brown	O	Transitions GEN 8 Emerald	E	Transitions XTRActive Polarised Gray	X
Remove 420 Photomatic Gray	U	Transitions XTRActive Brown	C	Photomatic Brown	B
Transitions GEN 8 Brown	Z	Transitions XTRActive Gray	L	Photomatic Gray	I
Transitions GEN 8 Gray	N	Transitions XTRActive Green	K	DriveWear	D
Transitions GEN 8 Green	V	Nupolar Brown	B		